

Associated with HISTO Stains Procedures Version Date: July 14, 2025

Page 1 of 1

|  | PDS#   |
|--|--|
| Histology Laboratory Special Stains Submissi                         | 1 Form - Contract and Research Cases                                     |
| Name of researcher:  |  |
| Contact name and phone #:  |  |
| Additional information:  |  |
| Note: If you have any questions, call the Histology La               | aboratory at 306-966-7319 or 306-966-7322.                               |
| Special Stain Request  |  |
| Aluminum - Aluminum  | Methyl Green Pyronine – plasma cells                                     |
| Alcian Blue – acidic mucins  | Oil Dad O limida Natau places submit wat tissue                          |
| Alcian Blue/PAS – acidic, neutral mucins<br>Alizarin Red S – calcium | Oil Red O – lipids – Note: please submit wet tissue Okajima – hemoglobin |
| Anzarin Red S – carcium<br>Congo Red – amyloid                       | Periodic Acid Methenamine Silver – basement                              |
| Congo Red – aniyloid<br>Duffy's – mast cells and eosinophils         | membranes  |
| Fites – Mycobacterium spp. and Nocardia spp.                         | Periodic Acid Schiff – neutral mucins, glycogen, fungu                   |
| Fontana Masson – melanin, late lipofuscins                           | Teriodic red beinir — neddal maems, grycogen, ranga                      |
|  | Periodic Acid Schiff's/Diastase – glycogen                               |
| Formalin Pigment Removal   | Perl's Prussian Blue – hemosiderin                                       |
| Giemsa – mast cell granules, bacteria                                | Phosphotungstic Acid Hematoxylin – fibrin, muscle                        |
| Gordon and Sweets – reticulin fibers                                 | striations   |
|  | Rhodanine – copper   |
|  | Safranine O/Fast Green – cartilage                                       |
|  | Schmorl's – lipofuscins, melanin   |
| Grocotts Methenamine Silver/H&E – fungus                             | Sirius Red (Picosirius) – osteoid  |
| Hall's – bile  | Toluidine Blue – mast cells  |
| H&E – routine stain  |  |
|  | Von Kossa – calcium salts  |
|  | Warthin Faulkner – spirochetes, bacteria                                 |
|  | Weigert Van Gieson – muscle, collagen, elastic fibers                    |
|  | Ziehl Neelsen – Mycobacterium spp.                                       |
|  | Ziehl Neelsen – Long – lead inclusion bodies                             |
| Luxol Fast Blue/H&E – myelin<br>Macchiavello's – Rickettsia spp.     | Other (please list)  |
| Masson's Trichrome – muscle, collagen fibers                         | <del></del>  |
| MILL DI LATOR 1  |  |
| Netami Bleach/HCE - Inclaim pigment temoval                          |  |
| For PDS use only   |  |
| PDS#: Date received:   | Date completed:  |
| Work relinquished to whom:(Name and                                  |  |
| Billing codes applied: Cost:   | Duie)  |
|  |  |

(Check all that apply) Blocks Slides Tissue forwarded to DSO (Date & Initials):\_

Billing done by (Date & Initials): \_

Charges checked by (Date & Initials): \_

## **End of Form**

Supersedes: August 22, 2022